

PULSE LAVAGE

2500K4-DP-C PULSATING SALINE LAVAGE, PRESSURIZED CO2 AND SUCTION





Pulsating saline lavage, pressurized Co2 and suction

Specifically made for orthopedic arthroplasty bone bed preparation.

The pulse lavage has the ability to significantly clean and remove debris, whilst simultaneously applying suction. Pressurized CO2 dries cancellous bone after the saline lavage, removing residual fluid and fatty material prior to cement implantation, and therefore diminishing excess fluid between the cement and bone. Increased bone cement penetration also increases implant fixation, thus reducing the arthroplasty revi-

Pulse Lavage Trigger «

Co2 Trigger

Battery Pack can be Disposable Type or Rechargeable Type



#K4-RB Rechargeable Battery Pack





#K4-CH

Battery Charger Pressurized CO2

Before Pressurized CO2

2500 K4-DPRB-C

sion rate.

Not included Battery Pack Rechargeable Battery Pack optional

2500K4-DP-C

Built-in Disposable Battery Pack

K4-RB

Rechargeable Battery Pack

Battery Charger (for K4-RB)

Advantages

- ✓ Built-in battery pack design without additional battery box
- ✓ Ergonomic design comfortable to handle and use
- ✓ Battery pack easy to remove user friendly design.
- ✓ Easy attached & released tip accessory
- ✓ Sterile contents, zero infection
- Tube length 3M \(\begin{align*}
 \)



& 886-4-2493-5023

886-4-2493-5032

info@kaiser-tech.com.tw

 info@kaiser-tech.com.tw

www.kaiser-tech.com.tw

• 4F, No. 1, Keji Road, Dali District, Taichung City 41264, TAIWAN.